

Bedside cabinet 2020

Manufactured from high quality materials

The frame of the cabinet table is white epoxy-coated steel. Front panel of the drawers and surface tops are solid birch or beech Iki high pressure laminate. The height of the meal tray is easy to adjust and it can be changed for left-handed and right-handed use without any tools.

The cabinet has four lockable wheels, a safety limiter to avoid the risk of the drawers being pulled out and a slipper rack as standard.



LOJER® For easy care

Technical data and features	
Width	49 cm
Depth	57 cm
Height	90 cm
Tray height adjustment range	60 - 100 cm
Weight	35 kg
Drawer dimensions	34 x 37 x 12 cm (2 pcs)
Meal tray	49 x 33 cm
Warranty	2 years
Country of Manufacture	Finland (Lojer Oy)

Factory installed options:

- Adapter for attaching nurse call device
- Lock for upper drawer

Other accessories:

- Dining tray that tilts to reading position
- Chromium plated side basket



LOJER® *For easy care*

© Lojer Group 2019. All Rights reserved | Lojer Group, Putajantie 42, 38210 Sastamala | Tel. + 358 (0) 10 830 6700 | myynti@lojer.com | www.lojer.com

Lojer Group is the largest manufacturer of medical furniture in the Nordic countries. In our two factories in Finland, we manufacture products such as hospital beds and operating tables, examination tables and treatment tables.